

BLK protein (STJ505671)
STJ505671**GENERAL INFORMATION**

Product Type	Proteins
Short Description	BLK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	640
Gene Symbol	BLK
Uniprot ID	BLK_HUMAN
Immunogen	Tyrosine-protein kinase BLK Protein
Immunogen Region	
Specificity	
Immunogen Sequence	