

Rabbit Anti-Goat Heavy Chain protein (STJ505674)

STJ505674

GENERAL INFORMATION

Product Type	Proteins
Short Description	Rabbit Anti-Goat Heavy Chain-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Rabbit anti-Goat IgG (H)
Immunogen Region	
Specificity	
Immunogen Sequence	