

B19 protein (STJ505678)
STJ505678**GENERAL INFORMATION**

Product Type	Proteins
Short Description	B19-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26559
Gene Symbol	Hunk
Uniprot ID	HUNK_MOUSE
Immunogen	Purified Hormonally Up-regulated Neu tumor-associated Kinase protein.
Immunogen Region	
Specificity	
Immunogen Sequence	