

Myosin 1C protein (STJ505679)

STJ505679

GENERAL INFORMATION

Product Type	Proteins
Short Description	Myosin 1C-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	17913
Gene Symbol	Myo1c
Uniprot ID	MYO1C_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1-50 on mouse Unconventional Myosin-Ic/Myo 1C.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	