

Myosin 1C protein (STJ505680)

STJ505680

GENERAL INFORMATION

Product Type	Proteins
Short Description	Myosin 1C-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	17913
Gene Symbol	Myo1c
Uniprot ID	MYO1C_MOUSE
Immunogen	Purified Unconventional Myosin-Ic/Myo 1C protein.
Immunogen Region	
Specificity	
Immunogen Sequence	