

RBP2 protein (STJ505683)

STJ505683

GENERAL INFORMATION

Product Type	Proteins
Short Description	RBP2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74152
Gene Symbol	Stra6l
Uniprot ID	STR6L_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse Retinol-Binding Protein 4 Receptor 2/RBP4 receptor 2/STRA6-Like Protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	