

TMC2 protein (STJ505687)

STJ505687

GENERAL INFORMATION

Product Type	Proteins
Short Description	TMC2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	117532
Gene Symbol	TMC2
Uniprot ID	TMC2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 850-906 on human Transmembrane channel-like protein 2.
Immunogen Region	850-906
Specificity	
Immunogen Sequence	