

CD39 protein (STJ505689)

STJ505689

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD39-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	953
Gene Symbol	ENTPD1
Uniprot ID	ENTP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 440-490 on human EctoNucleoside TriPhosphate Diphosphohydrolase 1.
Immunogen Region	440-490
Specificity	
Immunogen Sequence	