

CD39 protein (STJ505690)

STJ505690

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD39-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	953
Gene Symbol	ENTPD1
Uniprot ID	ENTP1_HUMAN
Immunogen	EctoNucleoside TriPhosphate Diphosphohydrolase 1 purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	