

Filamin B protein (Phospho) (2478) (STJ505691)

STJ505691

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-Filamin B-2478 protein was developed from and has a target region of 2478. For use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2317
Gene Symbol	FLNB
Uniprot ID	FLNB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 2450-2500 on human Phospho-Filamin-B. Phosphorylated 2478 serine residue.
Immunogen Region	2450-2500
Specificity	
Immunogen Sequence	