

DAP Kinase 1 protein (STJ505694)

STJ505694

GENERAL INFORMATION

Product Type	Proteins
Short Description	DAP Kinase 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1612
Gene Symbol	DAPK1
Uniprot ID	DAPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1370-1420 on human Death-associated protein kinase 1 protein.
Immunogen Region	1370-1420
Specificity	
Immunogen Sequence	