

## DAP Kinase 1 protein (STJ505695)

STJ505695

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	DAP Kinase 1-protein if for use in research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:2, 500-1:5, 000
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	1612
<b>Gene Symbol</b>	DAPK1
<b>Uniprot ID</b>	DAPK1_HUMAN
<b>Immunogen</b>	Purified Death-associated protein kinase 1 protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	