

GBA protein (STJ505702)

STJ505702

GENERAL INFORMATION

Product Type	Proteins
Short Description	GBA-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2629
Gene Symbol	GBA
Uniprot ID	GLCM_HUMAN
Immunogen	Recombinant Glucosylceramidase protein.
Immunogen Region	
Specificity	
Immunogen Sequence	