

Bile Acid Receptor NR1H4 protein (STJ505708)

STJ505708

GENERAL INFORMATION

Product Type	Proteins
Short Description	Bile Acid Receptor NR1H4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9971
Gene Symbol	NR1H4
Uniprot ID	NR1H4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on huma Bile acid receptor/NR1H4/FXR A3-4 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	