

Bile Acid Receptor NR1H4 protein (STJ505709)

STJ505709

GENERAL INFORMATION

Product Type	Proteins
Short Description	Bile Acid Receptor NR1H4-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9971
Gene Symbol	NR1H4
Uniprot ID	NR1H4_HUMAN
Immunogen	Bile acid receptor/NR1H4/FXR A3-4
Immunogen Region	
Specificity	
Immunogen Sequence	