

CD39 protein (STJ505711)

STJ505711

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD39-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12495
Gene Symbol	Entpd1
Uniprot ID	ENTP1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 380-430 on mouse EctoNucleoside TriPhosphate Diphosphohydrolase 1 protein.
Immunogen Region	380-430
Specificity	
Immunogen Sequence	