

P70 S6 Kinase beta protein (STJ505718)

STJ505718

GENERAL INFORMATION

Product Type	Proteins
Short Description	P70 S6 Kinase beta-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6198
Gene Symbol	RPS6KB1
Uniprot ID	KS6B1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Ribosomal protein S6 kinase beta protein.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	