

PICK1 protein (STJ505720)
STJ505720**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PICK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9463
Gene Symbol	PICK1
Uniprot ID	PICK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 330-380 on human PRKCA-binding protein/PICK1 protein.
Immunogen Region	330-380
Specificity	
Immunogen Sequence	