

IFI6 protein (STJ505722)
STJ505722**GENERAL INFORMATION**

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | IFI6-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 2537 |
| Gene Symbol | IFI6 |
| Uniprot ID | IFI6_HUMAN |
| Immunogen | Synthetic peptide take within amino acid region 80-130 on human Interferon alpha-inducible protein 6. |
| Immunogen Region | 80-130 |
| Specificity | |
| Immunogen Sequence | |