

SRPK1 protein (STJ505726)
STJ505726**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SRPK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6732
Gene Symbol	SRPK1
Uniprot ID	SRPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 405-455 on human SRSF protein kinase 1.
Immunogen Region	405-455
Specificity	
Immunogen Sequence	