

PP2A-alpha protein (STJ505728)

STJ505728

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | PP2A-alpha-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------------------|---|
| Gene ID | 5525 |
| Gene Symbol | PPP2R5A |
| Uniprot ID | 2A5A_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 410-486 on human Serine/Threonine-Protein Phosphatase 2A, B56-Alpha subunit. |
| Immunogen Region | 410-486 |
| Specificity | |
| Immunogen | |
| Sequence | |