

## PP2A-beta protein (STJ505730)

STJ505730

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	PP2A-beta-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	5526
<b>Gene Symbol</b>	PPP2R5B
<b>Uniprot ID</b>	2A5B_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 450-497 on human Serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Beta subunit.
<b>Immunogen Region</b>	450-497
<b>Specificity</b>	
<b>Immunogen Sequence</b>	