

PP2A-epsilon protein (STJ505734)

STJ505734

GENERAL INFORMATION

Product Type	Proteins
Short Description	PP2A-epsilon-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5529
Gene Symbol	PPP2R5E
Uniprot ID	2A5E_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-467 on human serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Epsilon subunit.
Immunogen Region	400-467
Specificity	
Immunogen Sequence	