

PP2A-epsilon protein (STJ505735)

STJ505735

GENERAL INFORMATION

Product Type	Proteins
Short Description	PP2A-epsilon-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5529
Gene Symbol	PPP2R5E
Uniprot ID	2A5E_HUMAN
Immunogen	Purified Serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Epsilon subunit protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	