

## PP2A-gamma protein (STJ505736)

STJ505736

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	PP2A-gamma-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5527</a>
<b>Gene Symbol</b>	<a href="#">PPP2R5C</a>
<b>Uniprot ID</b>	<a href="#">2A5G_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 390-439 on human Serine/Threonine-Protein Phosphatase 2A B56 Gamma1 subunit.
<b>Immunogen Region</b>	390-439
<b>Specificity</b>	
<b>Immunogen Sequence</b>	