

PP2A-gamma protein (STJ505737)

STJ505737

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | PP2A-gamma-protein if for use in research applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel. |
| Conjugation | |
| Purification | |
| Dilution | |
| Range | |
| Formulation | |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | 5527 |
| Gene Symbol | PPP2R5C |
| Uniprot ID | 2A5G_HUMAN |
| Immunogen | Purified Serine/Threonine-Protein Phosphatase 2A B56 Gamma1 subunit. |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |