

OR14I1 protein (STJ505744)

STJ505744

GENERAL INFORMATION

Product Type	Proteins
Short Description	OR14I1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	401994
Gene Symbol	OR14I1
Uniprot ID	O14I1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Olfactory Receptor, family 14, subfamily I, member 1 protein.
Immunogen Region	1-50
Specificity	
Immunogen	
Sequence	