

OR14I1 protein (STJ505745)

STJ505745

GENERAL INFORMATION

Product Type	Proteins
Short Description	OR14I1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	401994
Gene Symbol	OR14I1
Uniprot ID	O14I1_HUMAN
Immunogen	Purified Olfactory Receptor, family 14, subfamily I, member 1.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	