

CRALBP protein (STJ505750)

STJ505750

GENERAL INFORMATION

Product Type	Proteins
Short Description	CRALBP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6017
Gene Symbol	RLBP1
Uniprot ID	RLBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-317 on human Retinaldehyde-binding protein 1/CRALBP/RLBP1 protein.
Immunogen Region	260-317
Specificity	
Immunogen	
Sequence	