

OTOP1 protein (STJ505752)

STJ505752

GENERAL INFORMATION

Product Type	Proteins
Short Description	OTOP1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	133060
Gene Symbol	OTOP1
Uniprot ID	OTOP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 370-420 on human Otopetrin-1.
Immunogen Region	370-420
Specificity	
Immunogen Sequence	