

OTOP2 protein (STJ505754)

STJ505754

GENERAL INFORMATION

Product Type	Proteins
Short Description	OTOP2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	237987
Gene Symbol	Otop2
Uniprot ID	OTOP2_MOUSE
Immunogen	Otopetrin-2
Immunogen Region	
Specificity	
Immunogen Sequence	