

IL-34 protein (STJ505758)

STJ505758

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL-34-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146433
Gene Symbol	IL34
Uniprot ID	IL34_HUMAN
Immunogen	Interleukin-34
Immunogen Region	
Specificity	
Immunogen Sequence	