

JAK2 protein (STJ505762)

STJ505762

GENERAL INFORMATION

Product Type	Proteins
Short Description	JAK2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3717
Gene Symbol	JAK2
Uniprot ID	JAK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 980-1030 on human JAK2 protein.
Immunogen Region	980-1030
Specificity	
Immunogen Sequence	