

PDE3B protein (STJ505766)

STJ505766

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE3B-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5140
Gene Symbol	PDE3B
Uniprot ID	PDE3B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-310 on human PDE3B protein. Selective only for PDE3B variant.
Immunogen Region	260-310
Specificity	
Immunogen Sequence	