

ZAP70 protein (STJ505770)
STJ505770**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ZAP70-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7535
Gene Symbol	ZAP70
Uniprot ID	ZAP70_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-410 on human Tyrosine-protein kinase ZAP-70.
Immunogen Region	350-410
Specificity	
Immunogen Sequence	