

ASXL3 protein (STJ505772)
STJ505772**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ASXL3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	211961
Gene Symbol	Asxl3
Uniprot ID	ASXL3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 350-410 on human Tyrosine-protein kinase ZAP-70.
Immunogen Region	350-410
Specificity	
Immunogen Sequence	