

ASXL3 protein (STJ505773)

STJ505773

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | ASXL3-protein if for use in research applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel. |
| Conjugation | |
| Purification | |
| Dilution | |
| Range | |
| Formulation | |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|------------------|
| Gene ID | 211961 |
| Gene Symbol | Asxl3 |
| Uniprot ID | ASXL3_MOUSE |
| Immunogen | Purified ZAP-70. |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |