

## PBXIP1 protein (Phospho) (STJ505778)

STJ505778

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Phospho-PBXIP1-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">229534</a>
<b>Gene Symbol</b>	<a href="#">Pbxip1</a>
<b>Uniprot ID</b>	<a href="#">PBIP1_MOUSE</a>
<b>Immunogen</b>	Phosphorylated synthetic peptide taken within amino acid region 100-150 on mouse Pre-B-cell leukemia transcription factor-interacting protein 1.
<b>Immunogen Region</b>	100-150
<b>Specificity</b>	
<b>Immunogen Sequence</b>	