

PBXIP1 protein (Phospho) (STJ505779)

STJ505779

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PBXIP1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	229534
Gene Symbol	Pbxip1
Uniprot ID	PBIP1_MOUSE
Immunogen	Purified Phosphorylated Pre-B-cell leukemia transcription factor-interacting protein 1.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	