

ADA protein (STJ505780)

STJ505780

GENERAL INFORMATION

Product Type	Proteins
Short Description	ADA-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100
Gene Symbol	ADA
Uniprot ID	ADA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human Adenosine Deaminase/Adenosine Aminohydrolase.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	