

## ADA protein (STJ505781)

STJ505781

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	ADA-protein if for use in research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution</b>	
<b>Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	100
<b>Gene Symbol</b>	ADA
<b>Uniprot ID</b>	ADA_HUMAN
<b>Immunogen</b>	Purified Adenosine Deaminase/Adenosine Aminohydrolase.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	