

Serine/threonine-protein kinase 4/MST-1 protein (STJ505782)

STJ505782

GENERAL INFORMATION

Product Type	Proteins
Short Description	Serine/threonine-protein kinase 4/MST-1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6789
Gene Symbol	STK4
Uniprot ID	STK4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Mammalian STE20-like protein kinase 1/STE20-like kinase MST1/STK4.
Immunogen Region	400-450
Specificity	
Immunogen	
Sequence	