

STK3/MST-2 protein (STJ505785)

STJ505785

GENERAL INFORMATION

Product Type	Proteins
Short Description	STK3/MST-2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	WB: 1:500
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56274
Gene Symbol	Stk3
Uniprot ID	STK3_MOUSE
Immunogen	Purified Mammalian STE20-like protein kinase 2/STE20-like kinase MST2/STK3.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	