

MST3 protein (STJ505786)

STJ505786

GENERAL INFORMATION

Product Type	Proteins
Short Description	MST3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8428
Gene Symbol	STK24
Uniprot ID	STK24_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 460-517 on human Mammalian STE20-like protein kinase 3/STE20-like kinase MST3/STK24.
Immunogen Region	460-517
Specificity	
Immunogen Sequence	