

Hck protein (STJ505788)

STJ505788

GENERAL INFORMATION

Product Type	Proteins
Short Description	Hck-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3055
Gene Symbol	HCK
Uniprot ID	HCK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human Hematopoietic cell kinase/Tyrosine-protein kinase HCK.
Immunogen Region	350-400
Specificity	
Immunogen	
Sequence	