

TUFM protein (STJ505796)
STJ505796**GENERAL INFORMATION**

Product Type	Proteins
Short Description	TUFM-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7284
Gene Symbol	TUFM
Uniprot ID	EFTU_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 325-375 on human Elongation factor Tu, mitochondrial.
Immunogen Region	325-375
Specificity	
Immunogen Sequence	