

PI3R6 protein (STJ505798)

STJ505798

GENERAL INFORMATION

Product Type	Proteins
Short Description	PI3R6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146850
Gene Symbol	PIK3R6
Uniprot ID	PI3R6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 590-640 on human Phospholinoside 3-kinase regulatory subunit 6.
Immunogen Region	590-640
Specificity	
Immunogen Sequence	