

EHMT1/GLP1 protein (STJ505800)

STJ505800

GENERAL INFORMATION

Product Type	Proteins
Short Description	EHMT1/GLP1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79813
Gene Symbol	EHMT1
Uniprot ID	EHMT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 950-100 on human Histone-lysine N-methyltransferase EHMT1/GLP1.
Immunogen Region	950-100
Specificity	
Immunogen Sequence	