

Somatostatin Receptor 2 protein (STJ505803)

STJ505803

GENERAL INFORMATION

Product Type	Proteins
Short Description	Somatostatin Receptor 2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6752
Gene Symbol	SSTR2
Uniprot ID	SSR2_HUMAN
Immunogen	Purified protein.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	