

SARS-CoV-2 surface glycoprotein protein (STJ505808)

STJ505808

GENERAL INFORMATION

Product Type	Proteins
Short Description	SARS-CoV-2 surface glycoprotein-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken from unique amino acid sequence on surface glycoprotein [Severe acute respiratory syndrome coronavirus 2].
Immunogen Region	
Specificity	
Immunogen Sequence	